


# Pairing MemorialCare Omron EVOLV<sup>®</sup> – Blood Pressure Monitor

Only pair your MemorialCare blood pressure monitor after you've successfully connected your Virtual Exam Kit. You will only need to perform this setup once.

**Step 1:** Activate Bluetooth<sup>®</sup> on the mobile device or tablet where the TytoCare app is installed.

**Step 2:** Sign into the TytoCare app.


**Step 3:** Access the user menu from the top left of the TytoCare app and select **Omron BPM Pair**.



**Step 4:** On the next screen, follow the instructions to place the Omron EVOLV in pairing mode.



**Step 5:** Select **Pair** in the TytoCare app when prompted.



The Bluetooth<sup>®</sup> word mark and logos are registered trademarks owned by Bluetooth SIG, Inc.